SOL-MILLENNIUM®

Technical Data Sheet

Product specification



Troduct specification				
1. Product name	SOL-CARE™ Insulin Safety Syringe U-100 Insulin Only			
2. Description	The SOL-CARE™ 50 and 100 unit Insulin Safety Syringe U-100 Insulin Only is used to inject insulin into the body. The Insulin Safety Syringe U-100 Insulin Only is designed to aid in the prevention of needle stick injuries.			
3. Characteristics	After use, the health care professional fully depresses the plunger to engage the Luer Assembly. Once the Luer Assembly is engaged, pulling back the plunger causes the Adapter and the attached needle to be withdrawn into the safety of the barrel. In this position against the flange, lateral pressure on the plunger results in a controlled fracture of the plunger. Both the syringe and plunger are discarded in a Sharps container.			
	The Insulin Safety Syringe U-100 Insulin Only is used to inject insulin into the body. In			

addition, the Insulin Safety Syringe U-100 Insulin Only is designed to aid in the

5. Instructions for use

4. Intended use

See IFU

6. Sizes and REF numbers

REF	Size
100002B	0.5ml 29G*1/2"
100089B	0.5ml 30G*5/16"
100097B	0.5ml 30G*15/64"
100001IM	0.5ml 30G*1/2"
100002IM	0.5ml 29G*1/2"
100089IM	0.5ml 30G*5/16"
100095IM	0.5ml 30G*15/64"
100017B	1ml 29G*1/2"
100020B	1ml 30G*5/16"
100017IM	1ml 29G*1/2"
100081IM	1ml 30G*1/2"

prevention of needle stick injures.

REF	Size		
100020IM	1ml 30G*5/16"		
100071IM	1ml 28G*1/2"		
100094IM	1ml 30G*15/64"		
100091IM	1ml 31G*5/16"		

SOL-MILLENNIUM®

Technical information							
 List of Materials 		Component name		Material			
		Needle cap		HDPE:2911			
		Cannula			Stainless steel: SUS304		
		Needle hub			PP: 5250T		
		O ring			Medical Silicone		
		Locating ring			ABS: HI121		
		Gasket			TPE:181W-57		
		Plunger			PP: 5250T		
		Barrel			PP: 5250	Т	
		Locking ring			PP: 5250T		
		Adhesive			Loctite: 3321		
		Needle Lubricant			Silicon oil: Alpha 12500cst		
		Barrel Lubricant			Silicon oil: DC 360 12500cst		
2. Latex free		YES					
3. PHT / DEHP / PVC / BPA free		YES					
4. Shelf life		5 years					
5. Sterilization metho	od	Sterilized using Et	hylene Oxide				
6. Packaging	5.1 Sal	oc unit	0.5ml&1ml			100	Units per box
specification	o.i San	es unit	0.5ml &1ml			1000	Units per case
7. Technical Drawing	1	2 3 4 5	6 7 	*******	8	*******	9
	2. 3. 4. 5.	Needle cap Cannula Needle hub O ring Locating ring Gasket		7. Plu 8. Bar 9. Loc	_		

SOL-MILLENNIUM®

Quality and Regulatory information						
1. Quality certificate	Quality Management System according ISO 13485					
2. Product classification	Class IIa according to Annex IX of MDD 93/42/EEC					
	The product is compliant with the following standards and regulations:					
	Document reference	Title				
	ISO 8537:2016	Sterile single-use syringes, with or without needle, for insulin				
	ISO 7886-4:2006	Sterile hypodermic syringes for single use-Part 4: Syringes with re-use prevention feature				
	ISO 9626:1991/Amd1:2001	Stainless steel needle tubing for the manufacture of medical devices				
	ISO 10993- 1:2009/Cor1:2010	Biological evaluation of medical devices				
	ISO10993-4:2002	Biological evaluation of medical devices—part 4				
	ISO10993-5:2009	Biological evaluation of medical devices—part 5				
3. List of standards	ISO10993-7:2008/Cor1:2009	Biological evaluation of medical devices—part 7				
	ISO10993-10:2010	Biological evaluation of medical devices—part 10				
	ISO 15223-1:2012	Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied				
	ISO 11607- 1:2006/Amd1:2014	Packaging for terminally sterilized medical devices—Part 1				
	ISO 11607-2:2006	Packaging for terminally sterilized medical devices—Part 2				
	ISO 11737- 2:2006/Amd1:2014	Sterilization of medical devices-Microbiological methods- Part 2: Tests of sterility performed in the validation of a sterilization process				
	EN 980:2008	Symbols for use in the labelling of medical devices				
	EN 1041:2008+A1:2013	Information supplied by the manufacturer of medical devices				

REV	02	Date	04.10.2016