

Harris haematoxylin

REF. 361070

Dye for histo-cytology staining



IFU113A-RAL

For professional use only.
Please read all information carefully before using this device.

Table of contents

Intended use.....	1
Principle.....	1
Device description	2
Storage	2
Hazard classification and safety information.....	2
Personnel qualification	2
Specific equipment and reagents required but not provided	2
Operating procedure.....	3
Expected results.....	5
Performance.....	6
User quality Control.....	6
Other products.....	6
Recommendations, notes, and troubleshooting.....	6
Table of symbols and abbreviations	8
Bibliography.....	8
Change tracking	8

Intended use

Harris haematoxylin is intended to be used in combination with other staining devices for histo-cytology staining prior microscopic examination.

If applicable, RAL Diagnostics recommends using the associated RAL Diagnostics products and cannot guarantee that the expected results will be achieved if used in combination with products of other brands.

Principle

Harris Haematoxylin is a nuclear dye that usually associated to the cytoplasmic dyes OG6 and EA50 for Papanicolaou staining or with Shorr solution for Harris-Shorr staining.

Papanicolaou staining is the most used technique for the diagnosis of cervical cancer. In its regressive variation Harris Haematoxylin has a long-acting time (one to three minutes).

Papanicolaou staining is also perform for spermatozoa morphological study (spermocytogram).

Harris-Shorr staining is principally used in hormonal cytology and to carry out spermatozoa morphological study (spermocytogram).

Device description

Harris haematoxylin

Clear red violet solution


REF. 361070-0500	1 X 0.5 L
REF. 361070-1000	1 X 1.0 L
REF. 361070-2500	1 X 2.5 L

For a specific batch, refer to the analysis certificate of the batch available at my.ral-diagnostics.fr.

Storage

Storage temperature: 15-25°C away from light.

Bottle shelf life before and after opening: refer to the expiry date on the label.



Hazard classification and safety information

Harris haematoxylin

No labelling applicable

Personnel qualification

All samples and products must be handled by qualified and authorized personnel, using individual or collective protection, in accordance with the national directives in force in the laboratories. Personnel must also be aware of the classification of hazardous materials indicated on the label and the safety data sheet (available at my.ral-diagnostics.fr).

The specimen must be treated in accordance with procedures available in the laboratory and required by national authorities.

The diagnosis must be conducted by qualified and authorized personnel, in accordance with the procedures in force within the laboratory.

Specific equipment and reagents required but not provided

Absolute ethanol, sodium bicarbonate, magnesium sulfate, hydrochloric acid 37%, ammonia 20%, microscope slides and these following RAL Diagnostics devices:

OG6 Papanicolaou solution REF. 361630

EA50 Papanicolaou solution REF. 367600

Shorr Solution REF. 361100

This equipment may vary depending on the protocol. Please refer to the relevant protocol (see the section operating procedure) to ensure that you have the necessary equipment to carry out tests.

Operating procedure

The equipment used for sample processing must comply with the supplier's instructions for use.

Sample preparation

Specimen must treat in accordance with procedures available in the laboratory and promulgated by national authorities.

For spermocytogram make a smear et let it air dry.

If the smear was fixed by a cytofixer, remove cytofixer by dipping the slide in 50° ethanol during 20 to 30 minutes.

Reagents and instruments preparation

0,25%-hydrochloric acid solution: pour 6.757 mL of HCl 37% fill up to with 1 L of distilled water.

Ammonia alcohol solution: mix 3 mL of ammonia 20% and 97 mL of 70° ethanol

Acid-alcohol Solution: mix 300 ml of absolute ethanol with 2 ml of concentrated hydrochloric acid and 100 ml of distilled water.

Scott Solution: dissolve 3,5 g of sodium bicarbonate (NaHCO₃) and 20 g of magnesium sulfate (MgSO₄ 7H₂O) in 1 L of distilled water.

Protocols

The staining steps of the protocols indicated below consist of a successive dipping of the slides in the different staining baths.

Protocol for Papanicolaou regressive variation staining - Manual bath method - Manual microscopic analysis

Processing time: 17 min

Steps	Reagent	Time [mm: ss]	Indications
Hydrate	80° ethanol	00: 30	Can be extended to 1 min
Hydrate	70° ethanol	00: 30	
Hydrate	50° ethanol	00: 30	
Hydrate	Distilled water	00:30	
Stain	Harris haematoxylin	01:00	Can be extended to 3 min
Rinse	Distilled water	00:30	Can be extended to 1 min
Differentiate	0,25 hydrochloric acid solution	No	Dips six times
Rinse	Tap water	06: 00	No
Rinse	Distilled water	00: 30	No
Dehydrate	50° ethanol	00: 30	Can be extended to 1 min
Dehydrate	70° ethanol	00: 30	
Dehydrate	80° ethanol	00: 30	
Dehydrate	95° ethanol	00: 30	
Stain	OG6 Papanicolaou solution	01: 30	No
Rinse	95° ethanol	01: 00	2 x 30 sec. Can be extended to 2 min (2 x 1min)
Stain	EA50 Papanicolaou solution	01: 00	No
Rinse	95° ethanol	01: 30	3 x 30 sec. Can be extended to 3 min (3 x 1min)
Dehydrate	Toluene or xylene	No	Pass slide in
Mount	Toluene or xylene base mounting media	No	No

Protocol for spermocytograms Papanicolaou method - Manual microscopic analysis staining - Manual bath

Processing time: 35 min 03 s

Steps	Reagent	Time [mm: ss]	Indications
Fix	Ether-Alcohol	05: 00	Can be extended to 15 min
Dry	No	No	Open air
Hydrate	80° ethanol	00: 10	10 x 1 sec
Hydrate	70° ethanol	00: 10	
Hydrate	50° ethanol	00: 10	
Hydrate	Distilled water	00: 10	
Stain	Harris haematoxylin	03: 00	No
Rinse	Tap water	03: 00	Can be extended to 5 min
Differentiate	Acid-alcohol solution	00: 03	3 dips x 1 sec
Rinse	Tap water	03: 00	Can be extended to 5 min
Rinse	Scott Solution	04: 00	No
Rinse	Distilled water	00: 10	10 x 1 sec
Dehydrate	50° ethanol	00: 10	
Dehydrate	70° ethanol	00: 10	
Dehydrate	80° ethanol	00: 10	
Dehydrate	95° ethanol	00: 10	
Stain	OG6 Papanicolaou solution	02: 00	No
Rinse	95° ethanol	00: 20	2 series of 10 x 1 sec
Stain	EA50 Papanicolaou solution	05: 00	No
Rinse	95° ethanol	00: 10	2 series of 5 x 1 sec
Dehydrate	Absolute ethanol	02: 00	No

Dehydrate	Absolute ethanol / toluene 50/50 mix	01: 00	No
Dehydrate	Toluene or xylene	05: 00	Pass slide in
Mount	Mounting media	No	No

Protocol for Harris-Shorr regressive variation staining - Manual bath method - Manual microscopic analysis

Processing time: 03 min 36 s

Steps	Reagent	Time [mm: ss]	Indications
Hydrate	80° ethanol	00: 30	Can be extended to 1 min
Hydrate	70° ethanol	00: 30	
Hydrate	50° ethanol	00: 30	
Hydrate	Distilled water	00:30	
Stain	Harris haematoxylin	01:00	Can be extended to 3 min
Rinse	Distilled water	No	Pass in 2 baths
Differentiate	Ammonia alcohol solution	00: 06	6 dips x 1 sec
Rinse	Distilled water	No	Pass slide in
Dehydrate	70° ethanol	No	No
Dehydrate	95° ethanol	No	No
Stain	Shorr solution	00:30	Can be extended to 3 min
Rinse	95° ethanol	No	Pass slide in
Dehydrate	Absolute ethanol	No	Pass slide in
Dehydrate	Toluene or xylene	No	Pass slide in
Mount	Mounting media	No	No

Protocol spermocytograms by simplified Harris-Shorr staining - Manual bath method - Manual microscopic analysis

Processing time: 18 min 42 s

Steps	Reagent	Time [mm: ss]	Indications
Fix	70° ethanol	05: 00	No
Rinse	Tap water	00: 12	12 x 1 sec. Can be extended to 15 sec (15 x 1 sec)
Stain	Harris haematoxylin	03: 00	Can be extended to 3 min
Rinse	Tap water	00: 12	12 x 1 sec. Can be extended to 15 sec (15 x 1 sec)
Differentiate	Ammonia alcohol solution	00: 06	6 dips x 1 sec
Rinse	Tap water	00: 12	12 x 1 sec. Can be extended to 15 sec (15 x 1 sec)
Dehydrate	70° ethanol	01: 00	No
Dehydrate	95° ethanol	01: 00	No
Stain	Shorr solution	05: 00	No
Rinse	95° ethanol	01: 00	Pass slide in
Dehydrate	Absolute ethanol	02: 00	Get rid of the excess stain rubbing the slide with a compress
Dry	No	No	Pass slide in

Expected results

Papanicolaou regressive variation staining

Nuclei: more or less dark violet-blue

Eosinophilic cells cytoplasm: pink, sometimes red-pink or orange

Cyanophilic cells cytoplasm: blue, sometimes greenish

Orangeophilic cells cytoplasm: brilliant orangey

Spermocytograms Papanicolaou staining

Head piece -Nucleus: purple

Head piece-Acrosome: green blue

Flagellum: green

Midpiece: pale green

Assess in percentage:

- abnormalities of the head, midpiece and flagellum
- agglutinates
- leukocytes, erythrocytes, cells

Harris-Shorr regressive variation staining

Nuclei: more or less dark violet-blue

Eosinophilic cells cytoplasm: orangey to brown

Cyanophilic cells cytoplasm: pale green to dark green

Spermocytograms by simplified Harris-Shorr staining

Head piece -Nucleus: dark violet

Head piece-Acrosome: light violet

Flagellum: greyish blue

Midpiece: violet

Cytoplasmic Remains: red

Assess in percentage:

- abnormalities of the head, midpiece and flagellum
- agglutinates
- leukocytes, erythrocytes, cells

If observed results vary from those expected, please contact RAL Diagnostics technical service through your usual supplier for assistance.

Performance

This medical device is state of the art. Its analytical performance, scientific validity and medical relevance are assessed in the CE marking review.

To ensure product performance, use clean and dry laboratory equipment.

The laboratory is responsible for notifying the manufacturer and state competent authority of any serious incident relating to the use of the medical device.

User quality Control

Users are responsible for determining the appropriate quality control procedures for their laboratory and complying with applicable laboratory regulations.

RAL Diagnostics recommend quality control at reagents renewal and for the first staining cycle of each day. Slides stained for quality control purposes should be checked to ensure that they are satisfactory for intended test (properly stained and free of precipitate). Staining results for each cell type must also be compliant with this manual expected results.

These quality control procedures should only be performed by qualified personnel.

Other products

For more information contact your usual supplier.

Recommendations, notes, and troubleshooting

Products appearance

If the appearance of the products differs from the description above, do not use it and contact RAL Diagnostics technical service through your usual supplier for assistance.

Procedures notes

To prevent products degradation, please comply with the storage and handling recommendations specified in this manual.

Reagent may present crystals. RAL Diagnostics recommend filtering the necessary amount of product for the staining before use.

Staining times may vary according to the nature of the smear or may be modified depending on the frequency of use, the desired staining intensity and the staining material used.

To avoid any detachment of the smear for the spermocytograms application, make a smear that would be neither too fine, nor too thick and allow it to dry well (during several hours in air or on a hot plate).

According to the thickness of the smear, it may be necessary to increase the fixative time in a bath of 70° ethanol, to 1 hour.

Products stability

Every RAL Diagnostics product can be used until the expiry date indicated on, in its original packaging if it is still hermetically sealed.

Staining stability

Staining quality and reproducibility depend on the correct use of the products.

RAL Diagnostics recommends mounting the stained slides with a coverslip using a suitable mounting liquid and to store them in a light and dustproof container.

Instructions for cleaning and waste disposal

All biological samples, effluents and used consumables should be considered potentially hazardous.



To avoid any risk, apply the following instructions: dispose of samples, effluents and consumables in accordance with laboratory standards and applicable national and local standards and regulations.

Chemical and biological waste must be collected and processed by specialized, registered companies.

Table of symbols and abbreviations

Depending on the product, you may find the following symbols on the device or the packaging material.

GHS PICTOGRAMS	INTERPRETATION	SYMBOL	INTERPRETATION
	Explosive		Batch code
	Flammable		Serial number
	Oxidizer		Catalogue reference
	Compressed gas		Date of manufacture
	Corrosive		Use up to
	Toxic		Unique device identifier
	Harmful		Manufacturer
	Health Hazard		Importer
	Environmental Hazard		Entry distributing the medical advice in the region concerned
	No labelling applicable		CE marking device
			In vitro diagnostic medical device
			Authorised Representative in the European Community
			Authorised Representative in Switzerland
			Complies with UK guidelines
			Do not use if packaging is damaged
			Keep away from light
			Temperature limit: 15-25°C
			Temperature limit: 15-30°C
			Keep dry
			Box: handling upwards
			Fragile
			Sterilised by irradiation
			Single sterile barrier system with outer protective packaging
			Sterile and radiation-sterilised barrier suit
			Do not reuse
			Do not resterilize
			Contents sufficient for n tests
			Hazardous material contained
			Consult instructions for use
			Use
			After opening, use within XX months
			The product must not be used in conjunction with an automatic colouring machine
			Indicates a medical device that contains potentially carcinogenic, mutagenic or reprotoxic (CMR) substances, or substances classified as endocrine disruptors

Bibliography

BOON M.E., DRIJVER J.S., *Routine Cytological Staining Techniques*, éd. MACMILLAN EDUCATION LTD, 1986, p. 71-79 et 129-130.

CASTELAIN G., *Technique du Cytodiagnostic*, Réactifs RAL, p. 5-7.

DUBOST R., *Technique du spermocytogramme*, Pharm. Biol., vol. 13, 1979, p. 133-134.

GANTER P., JOLLES J., *Histochimie normale et pathologique*, éd. Gauthier-Villars, 1970, vol. 2, p. 1477.

LIENARD G., *Techniques d'investigation du sperme, spermogramme et spermocytogramme* (service du Professeur Salesse, Hôpital Cochin), UNATEB, « Thème et Débat » du 11 mars 1980.

HOULD R., *Techniques d'histopathologie et de cytopathologie*, éd. Maloine - Décarie, 1984, p. 312-313.

WHO Laboratory Manual, *For the examination of human semen and sperm-cervical mucus interaction*, 4^{ème} éd., World Health Organization, 1999, p. 75.

Procedure modified by Service de Biologie de la Reproduction, Hôpital Pellegrin - Bordeaux - FRANCE.

Change tracking

Date	Version	Changes
05/2022	IFU113A-RAL	IVDR (EU) 2017/746 compliance

