Technical sheet



Safranin

Code product : **361500**



Appearance: Absorption Maximum in water:

Clear red solution 515 nm $\leq \lambda \leq$ 520 nm

Applications:

Bacteriology: used in Gram-Hücker staining for the detection of Gram-positive and Gram-negative bacteria. Safranin stains Gram-negative bacteria pink.

Parasitology and Mycology: used in Kit Gram-Hücker L (slow action) for the detection of pathogenic fungi as well as microsporidia spores (protozoa).

Updated: 25/04/2013