





# Combination of epT.I.P.S.<sup>®</sup> pipette tips with Xplorer/Xplorer plus and Xplorer plus Move It<sup>®</sup> pipettes

Eppendorf Xplorer <sup>®</sup> /Eppendorf Xplorer <sup>®</sup> plus/Eppendorf Xplorer <sup>®</sup> plus Move It <sup>®</sup>	0.1 – 10 µL, dark gray	0.1 – 10 µL, M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL, L, light gray	0.1 – 20 µL, light pink	0.5 – 100 µL, light yellow	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow
epT.I.P.S. <sup>®</sup>	■		■	■	■	■	■	■	■
epT.I.P.S. <sup>®</sup> LoRetention	■			■					■
ep Dualfilter T.I.P.S. <sup>®</sup>	■	■		■	■	■	■	■	■
ep Dualfilter T.I.P.S. <sup>®</sup> LoRetention	■			■				■	
ep Dualfilter T.I.P.S. <sup>®</sup> SealMax	■			■				■	■
<b>Variable volume, single-, 8- and 12-channel pipettes</b>									
0.5 – 10 µL	■	■	■	■					
1 – 20 µL	■	◆	◆	■					
5 – 100 µL	■						◆	■	■
10 – 200 µL	■						◆	◆	■
15 – 300 µL	■						◆	◆	◆
50 – 1,000 µL	■								
50 – 1,200 µL	■								
0.1 – 2.5 mL	■								
0.2 – 5 mL	■								
0.5 – 10 mL	■								
<b>Variable volume, 16- and 24-channel pipettes</b>									
1 – 20 µL	■				■				
1 – 20 µL	■					■			
<b>Move It<sup>®</sup>, variable volume, 4, 6, 8 and 12-channel pipettes</b>									
1 – 20 µL	■				■				
5 – 100 µL	■					■			
15 – 300 µL	■						◆	◆	◆
50 – 1,200 µL	■								

■ normal volume ◆ limited volume

20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL, L, dark green	0.25 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise	GEloader <sup>®</sup> /Microloader™ 0.5 – 20 µL
■	■	■	■	■	■	■	■	■	■
■	■								
■	■								
■	■								
■	■								
									■
◆									
■									
■									
	■	■	■						
		■	■						
				■					
					■	■			
							■	■	
■									
		■	■						