

## Crystal Violet Oxalate

Code product: **361490**

### Specifications:

- **Aspect:**
- **Absorption Maximum in water:**

Clear dark violet Solution  
 $588 \text{ nm} \leq \lambda \leq 592 \text{ nm}$

### Applications:

**Bacteriology:** used in Gram-Hücker staining for the detection of Gram-positive and Gram-negative bacteria. Crystal Violet Oxalate stains Gram-positive bacteria violet.

**Parasitology and Mycology:** used in Kit Gram-Hücker staining with Slow Differentiator for the detection of pathogenic fungi as well as microsporidia spores (protozoa).

**Histo-Cytology:** used in Gram staining according to Brown and Brenn for the detection of Gram-positive and Gram-negative germs on histological sections.

Updated: 25/04/2013