



Biosigma SpA

a Dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.)
| Fax ++39 0426 302228 | SMS ++39 348 4077376 |
E-mail info@biosigmaeu.com | <http://www.biosigma.com> |

M521D

Ultimo aggiornamento
30/05/2023

SCHEDA TECNICA PRODOTTO/ TECHNICAL DATA SHEET

REF: BSA0594

Provetta 13x75 mm, 6 ml senza bordo in PS / Test tube 13x75 mm, 6 ml without rim in PS

REF: BSA0594/P

Provetta 13x75 mm, 6 ml senza bordo in PP / Test tube 13x75 mm, 6 ml without rim in PP



BSA0594 Polistirolo/ Polystyrene



BSA0594/P Polipropilene/ Polypropylene

PER USI GENERICI DI LABORATORIO / FOR GENERAL LABORATORY USES

Famiglia: PROVETTE e TAPPI Family: TEST-TUBES & STOPPERS						
REF	Descrizione Description	Dimensioni Dimension	Conf. Box	Autoclavabile Autoclavable	Sterile	Vel. centrifug. Speed centrifugation
BSA0594	Provetta 6 ml in PS neutra Test Tube 6 ml cylindrical in PS neutral	Ø 13 x 75mm	4 x1.000	No/No	No/No	4.000 rpm
BSA0594/P	Provetta 6 ml in PP neutra Test Tube 6 ml cylindrical in PP neutral	Ø 13 x 75mm	4x1.000	Si/Yes	No/No	4.000 rpm

Destinato esclusivamente ad uso professionale

RACCOMANDAZIONI D'USO (a carattere esclusivamente informativo)

- In caso di fuoriuscita del contenuto usare i guanti pericolo di contaminazione
- Non utilizzare il prodotto scaduto
- Non variare la destinazione d'uso
- Prodotto non adatto ai bambini
- Conservare in luogo asciutto, lontano da fonti di calore e dalla luce solare diretta
- Smaltimento: utilizzare gli appositi D.P.I. e smaltire secondo la normativa vigente

For professional use only

RECOMMENDATION USE (exclusive as information)

- Use gloves: contamination risk if contents leaks
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place and far from sources of heat and from direct sunlight
- Waste disposal: use the provided D.P.I. and discharge under the laws in force.

Definizione del fabbricante / Manufacturer definition



BIOSIGMA SpA

Via Valletta, 6 - 30010 - Cantarana di Cona (VE) - ITALY
Tel. +39.0426.302224 - Fax +39.0426.302228
E-mail: info@biosigma.com - <http://www.biosigma.com>