

	Compatibiliteitstabel ART-tips				
Model	Gilson	Finnpipette	Biohit	Eppendorf	Socorex
ART 10 Reach	P2-P10	2 en 10 µl micro	10 µl	2.5 µl	
ART 10	P2-P10	2 en 10 µl	10 µl		
ART 20E				10 µl	10 µl
ART 10F		10 µl universeel	10 µl		
ART 20P	P20			20 µl	20 µl
ART 20				10 µl	
ART 50U		5 - 50 µl (multikanaal)	5 - 50 µl (multikanaal)	5 - 50 µl (multikanaal)	
ART 100		40 µl, 5 - 50 µl (multikanaal)	5 - 50 µl (multikanaal)		
ART 100E	P100		50 µl en Electronic 100 µl	100 µl	
ART 200	P200		200 µl	200 µl	
ART 200U	P200	50 - 300 µl (multikanaal), 200 µl	200 µl	50 - 300 µl (multikanaal), 200 µl	200 µl
ART 200G	P200	200 µl	200 µl	200 µl	
ART 300		300 µl (multikanaal)	250 µl (multikanaal)	Research Pro 300 µl	
ART 500			50 - 500 µl	50 - 250 µl	
ART 1000	P1000	1000 µl	1000 µl	1000 µl	
ART 1000E	P1000	1000 µl	1000 µl		
ART 1000G	P1000	1000 µl	1000 µl	1000 µl	
ART 1000 Reach	P1000	1000 µl	1000 µl	1000 µl	
ART 5000	P5000				
ART GEL 20P	P20	10 µl		20 µl	
ART GEL	P100	40 µl		100 µl	
ART XLP	P100	40 µl		100 µl	
ART XLP 200	P200			200 µl	1000 µl
ART XLG	P200	40 µl		100 en 200 µl	