

Technical Sheet

EOSIN-RAL 555

Product code: 361640

Specifications:

- **Appearance:** Clear orangey red Solution
- **Absorption Maximum in Water:** $516 \text{ nm} \leq \lambda \leq 518 \text{ nm}$
- **pH:** 6.9 – 7.1

Applications:

Haematology: used in Kit RAL 555, fast variation of May-Grünwald Giemsa staining, for differential staining of medullary and blood smears.

Parasitology and Mycology: used in Kit RAL 555 to detect parasites and fungi on blood smears.

Histo-Cytology: used in Kit RAL 555 for fast polychromic staining of histological sections and for the detection of *Helicobacter pylori* on histological sections.

Updated: 07/02/2022