

























## Socorex compatibility chart

Chart indicates optimal combination of Qualitix® and Socorex pipettes.

		Acura® manual and Acura® electro																						
		815					825 / 826 / 926					835 / 936			855 / 956									
		Up to 10 µL	10 y/15/20 µL	25 to 100 µL	120 to 200 µL	250 to 1000 µL	01 - 2 µL	0.5 - 10 µL	1 - 10 µL (*)	2 - 20 µL	5 - 50 µL	10 - 100 µL	20 - 200 µL	100 - 1000 µL	0.2 - 2 mL	0.5 - 5 mL	1-10 mL	0.5 - 10 µL	2 - 20 µL	5 - 50 µL	10 - 100 µL	20 - 200 µL	40 - 350 µL	
<b>10-20 µL ultra and microtips</b>	<b>Tip ref.</b>																							
	3090	●					●	●										●						
	3020	●					●	●											●					
	302	● <sup>2</sup>						● <sup>2</sup>																
	3090	●					●	●										●						
	3020	●					●	●											●					
	307	●					●	●										●						
<b>20 - 100 µL Tips</b>	<b>Tip ref.</b>																							
	3080		●					●	●	● <sup>1</sup>										● <sup>1</sup>				
	3080			●						●	●	● <sup>1</sup>								●	●	● <sup>1</sup>		
<b>200 µL Tips</b>	<b>Tip ref.</b>																							
	3070		●	●	●			●	●	●	●	●								●	●	●		
	3290		●	●	●			●	●	●	●	●								●	●	●		
	3080				●				●	●	●	●												
	3070			●	●							●											●	
	3080				●							●											●	
	307/327		●	●	●			●	●	●	●	●								●	●	●		
<b>300 - 350 µL Tips</b>	<b>Tip ref.</b>																							
	308				●						●	●												●
	308				●							●												● <sup>1</sup>
<b>1000 µL Tips</b>	<b>Tip ref.</b>																							
	3070				●								●											
	3190				●								●											
	3090				●								●											
	3070				●								●											
	307/317				●								●											
<b>2, 5, 10 mL Macrotips</b>	<b>Tip ref.</b>																							
	312													●										
	312														●									
	312															●								

<sup>1</sup> Up to max. tip volume    <sup>2</sup> For volumes > 2 µL

