



Kit Trichrome de Masson, Aniline Blue Variation

Code product : **361670-0000**

Composition:

• Mayer Haemalum	1 x 100 ml
• Ponceau Fuchsin	1 x 100 ml
• 1% Acetic Water	2 x 100 ml
• Phosphomolybdic Acid	1 x 100 ml
• Aniline Blue	1 x 100 ml

Applications:

Histo-Cytology:

Staining with kit Trichrome de Masson, Aniline Blue Variation, enables a trichromic staining on histological sections with a specific coloration of the collagen thanks to Aniline Blue.

Mayer Haemalum stains nuclei black-blue or brown.

Ponceau Fuchsin Solution stains cytoplasms red-pink.

Aniline Blue stains collagen blue.

Updated: 25/04/2013