

Pluggo Abbott Alinity 6 Position Automated Decapper - Technical Description (Draft)

Pluggo is designed for the safe removal of original caps from most standard vacuum collection tubes. Pluggo can handle 5ml, 7ml and 10ml vacuum tubes with 75 x 13 mm, 100 x 13 mm, 100 x 16 mm tubes carried in the Alinity racks.



PLUGGO

A detachable chute guides the caps into a disposal area after decapping.

Techical data

Capacity: Sample tubes: Speed: Size: Weight: Voltage: Power: Up to 6 samples in each rack Standard blood collection tubes 8 seconds per rack 29 x 27 x 21cm (Depth x Width x Height) 9,0 kg 230VAC 50Hz / 110VAC 60Hz 60VA

🖂 info@clinicon.co.uk	L +44 (0) 1273 463609	www.clinicon.co.uk