



YVSOLAB NV
Veedijk 33
2300 TURNHOUT
BELGIUM

Tel:+32 (0)14/67 27 79
Fax:+32 (0)14/55 77 68
E-mail: info@yvsolab.be
Website: www.yvsolab.be

TECHNICAL DATA SHEET

F.A.A. fixative (AFA)	
Product code: SL.AFAP.XX/XX // RE.2BAXX.XX/XX Date of issue: 23/11/2017 (v5)	
Expected aim	Product for the preparation of cyto-histological samples for optical microscopy.
Application	Fixative mixture used for samples with high fat tissue content.
Components	Acetic acid; Ethanol Denat;. Demi-water;
Warning and precaution	The product must be used exclusively by specialized technical operators. Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risks presented by the mixture, the precautionary measures during use, the measures first aid and the intervention in the event of accidental release. Do not use if the primary container is damaged. The safety data sheet is available on request.
Storage	Store the preparation at room temperature. Keep the containers tightly closed.
Stability	After the first opening, the product is usable until the expiry date, if correctly stored.
Disposal	Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.
Information in this publication is believed to be accurate and is given in good faith, but it is for the customer to satisfy itself of the suitability for its own particular purpose. No representation, warranty or guarantee is made as to its accuracy, reliability or completeness.	
 	