



Kit RAL 555

Code product: 361550-0000

Composition:

• FIX-RAL 555	Code 362870-0100	1 x 100 ml
• EOSIN-RAL 555	Code 361640-0100	1 x 100 ml
• BLUE-RAL 555	Code 361650-0100	1 x 100 ml

Applications:

Haematology:

Kit RAL 555, fast variation of May-Grünwald Giemsa technique, enables a differential staining of medullary and blood smears.

Cyto-Bacteriology:

Kit RAL 555 allows cytological study of urines, liquids and spinal fluid.

Parasitology and Mycology:

Kit RAL 555 enables to detect parasites and fungi on blood smears: tissular protozoa (*Leishmania*, *Toxoplasma*, Microsporidia), *Cryptosporidium*, *Pneumocystis carinii*, fungi contributing to deep mycosis, veterinary parasitology (Piroplasms, *M. pachydermatis*), *Plasmodium*, *Trichomonas*, Microfilaria.

Histo-Cytology:

Kit RAL 555 allows cytological study of fixed and paraffin embedded tissue sections, punctions (breast and deep organs), effusion fluids of serous membranes (pleura, peritoneum, ...), and enables to detect *Helicobacter pylori* on histological sections.

Updated: 25/04/2013