

TECHNICAL FEATURES

Trade Name	Graduated narrow neck bottles – without caps		
Use	An ideal bottle for water and oils		
Material	Polyhethylene HDPE		
Other features	There are no shoulders to restrict the flow		

DIMENSIONAL FEATURES

Heigh	mm. 271,00		
Base diameter	mm. 120,00		
Mouth diameter	mm. 23 (DIN GL 32)		
General dimesional tolerance	mm. +/- 1,00		
Wall tickness (bottle)	mm. 0,60 +/- 0,2		
Capacity	ml. 2000		



GENERAL FEATURES

CE Mark	No		
Foodstuff suitability	Yes		
Neck and thread	Conformity to DIN 13316 and 168		
Graduation	Graduated. Can be used as rough measures		
Sealing	Can be secured using shoulder and cap tag to avoid sample		
	contamination in transiti		

STERILIZATION OF THE BOTTLE

Autoclavable at +121°C	No
Gas (Ethylene Oxide)	Yes
Dry at +160°C	No
Chemical (Formalin)	Yes
Radiaions	Yes
Microwaves	Yes

Quality Assurance Massimiliano Capitanio



CLEANING OPERATIONS

Slight contamination	Wash with chemically neutral (pH7) cleaning agents
Heavy contamination	Wash with alkaline (pH up to 12) cleaning agent

CHEMICAL RESISTANCES

Categories of substances	Reference	Concent. %	T. 20° C	T. 40° C	T. 60° C
Inorganic Acid	Sulphuric acid	98	В	В	В
Organic Acid	Benzoic acid	100	А	А	А
Alcohol	Ethanol	100	А	В	В
Aldehydes	Acetaldeyde	100	В	С	-
Inorganic alkali	Isodium hydroxide	50	С	-	-
Organic alkali	Aniline	100	А	В	С
Ketone	Acetone	100	В	С	-
Esther	Ethyl acetate	100	В	С	-
Hydrocarbons, halogenated	Dichloroethylene	100	С	-	-
Hydrocarbons, aromatic	Benzene	100	В	С	-
Hydrocarbons, linear	Hexane	100	С	-	-
Supersolvent	Tetrahydrofuran	100	В	В	С

Legenda (not applicable if neglected):

- A: fair resistence; exposure (30 days) to the chemical does not cause any damage.

- **B: sufficient resistence;** exposure causes damage of poor importance, which sometimes is only temporary.

- C: poor resistence; exposure to chemical is not allowed, or causes immediate permanent damages.

Results of testing with reference materials are to be considered as indication: in case of specific use it is recommended to carry out preliminary testing.

Information about resistance of plastic materials to temperatures, steriization and washing treatments are based on bulletins issued by the Producers of raw materials , on literature data and on the experience gained in using the products

PACKAGING

Bag

15 pcs

Quality Assurance Massimiliano Capitanio Technical Direction Guido Borona