




LP ITALIANA SPA

Via C. Reale, 15/4 - 20157 Milano - ITALY
 Tel. +39 02 3933061
 Fax +39 02 39313484
 www.lpitaliana.com
 info@lpitaliana.com

Capitale Sociale € 309.600,00
 R.E.A. MI 882798
 Reg. Imp. MI 161285/3927/35
 C.F. e P.I. 01794050151
 c/c postale 19643204


TECHNICAL DATA SHEET 15 ml TEST TUBES

15 ml LP test tubes, with conical bottom and screw cap, are manufactured in biologically inert clear polypropylene.



Features / technical data:

- suitable for the collection, transport and test of every type of biologic sample (suspension cell, urine, bronchial excretion, etc).
- indicative graduation every 0,5 ml to 14 ml.
- materials: tube is made of virgin, semi-transparent polypropylene, with excellent chemical stability; items 111548 - 118548 of gamma-resistant PP; both items colour: natural; screw cap is made of high density polyethylene; colour: blue.
- the drawing and the thickness of the walls allows safe centrifugation up to 4.600 xg.
- tube wide temperature working range: from -5°C up to +121°C (freezing and autoclave).
- the leak-proof screw cap guards the operator from sample overflow and has temperature working range from -18°C and +80°C.
- items 111547 - 111548 are packed in bags of 125 pcs.
- item 118548 is on tray of expanded PS - 50 pieces each
- in the sterile version the sterility is assured until the bag is opened; internal sterility is guaranteed until the cap is removed.



LP codes and descriptions:

Code	111547	111548	118548	119548
Nominal Capacity ml	15			
Nominal Dimensions mm	17 x 120			
Sterility R	No	Yes		
Quantity per sealed bag	125	50	1	
Quantity per box	1000	200	250	

Intended use and CE mark: LP 15 ml test tubes intended use is to transport and keep fluids derived from the human body for the purpose of providing information concerning a physiological or pathological state of one patient, furthermore, where applicable, having a special microbiological state (indicated by **STERILE R** mark) they fall under the European Directive 98/79/CE - Annex III - on In Vitro Medical Diagnostic Devices (IVD), according to which they bear the conformity **CE** mark.

Disposal: LP 15 ml test tubes, screw cap, are totally recyclable; can also be incinerated while respecting antipollution norms because, after burning, they liberate, in the environment, water and carbon dioxide only (combustion performed at right temperature and with correct air injection).

Sterility: items 111548 - 118548 are sterilized by irradiation with ionizing radiations in conformity of one validated method, according to UNI EN ISO 11137-1:2006 and UNI EN ISO 11137-2:2012 at 10⁻⁶ SAL (Sterility Assurance Level), assured till opening of the packaging.

Five-year "shelf life", packed in sealed bag.

Firma:	Distribuito a: DG - DVE - VEE - VEI - UA		
Documento: Rssg 088 Rev. G	Fg. 1 di 1	Data Emissione: 01.07.2014	