



Export and Production
- since 1949 -

Heinz Herenz Medizinalbedarf GmbH • Rudorffweg 10 • 21031 Hamburg, Germany

Product description

Item:	Steristoppers [®] , cellulose stoppers
Material:	2-ply cellulose made of 100% elementary-chlorine-free cellulose
Sinking time:	max. 5 sec. (DAB 2010)
Ability of water holding:	min. 13 g/g (DAB 2010)
Glue:	Glue on basis of methyl-cellulose sprayed on the cellulose while twisting it
Shape:	Mould twisting process
Frequency of use:	for single use
Measuring tolerances:	due to softness of the material
Sterilization:	autoclavable or dry heat up to 180-200°C
Aspiration:	Steristoppers [®] are air permeable
Experience:	Steristoppers [®] have been used to seal tissue cultures for many years in bacteriological institutes all over the world to entire satisfaction
Patents:	German Patent No. 1467937 and foreign patents registered.
Manufacturing Lot No.:	printed on each carton: week / year
Shelf Life:	N/A

www.herenz.de

Phone +49-040-739-20-40
Fax +49-040-739-20-439
Email hhh@herenz.de
ILN 40 49434 00000 8
Steristopfen-e

Managing Directors

Dieter Cyll, Antje Herenz-Cyll
Comm.Reg. Hamburg HRB 48884
Tax No. 44/780/02888
V.A.T. No. DE 118 665 646

Postbank Hamburg

Sort Code 200 100 20
Account No. 6846 46 200
IBAN DE97200100200684646200
Swift-Code PBNKDEFF

Deutsche Bank AG

Sort Code 200 700 00
Account No. 55 00 251
IBAN DE82200700000550025100
Swift-Code DEUTDEHH

Hamburger Sparkasse

Sort Code 200 505 50
Account No. 1085 243 002
IBAN DE87200505501085243002
Swift-Code HASPDEHH

01.11.2017 rev. 2